

GM Series -Portable Balance
0.01g,0.1g,1g



Features	Benefit
Entry level model with lab precision	Cost effective and high performance
Ergonomically optimized key pad w/ large keys	Convenient and efficient performance
Precise tare button	Calculates the net weight of your ingredients
Secure and non-slip positioning with rubber feet	Precise leveling of the scale
Detachable pan	For easy and hygienic cleaning
Both AC and rechargeable batteries	Complete mobility
Second display or Rear display (Optional)	Convenient for customer to read the w

Model	GMB 302	GMB 602	GMB 1002	GMB 2002	GMB 3002	GMC 1002	GMC 2002	GMC 3002	GMC 5002	GMD 3001	GMD 6001	GME 6000
Draft Shield	With Draft Shield					No Draft Shield						
Capacity	300g	600g	1000g	2000g	3000g	1000g	2000g	3000g	5000g	3000g	6000g	6000g
Resolution	0.01 g									0.1 g		1 g
Repeatability	± 0.02 g								± 0.03 g	± 0.2 g		± 2 g
Linearity Error	± 0.02 g								± 0.03 g	± 0.2 g		± 2 g
Response Time	≤ 2 sec											
Calibration	Calibration with External Weight											
Weighing Units	g, ct, lb, oz, ozt, dwt, kg, PCS, GSM, TMR, PKT, etc											
Construction	static resistant ABS, stainless steel pan											
Display	Large, easy to read LCD display with backlight											
Pan Size	Round Pan (150mm)					Round Pan (180mm)						
Body dim.	230x 270x 70 mm (L x W x H)											
Battery	Rechargeable battery 6V 2.5A or Lithium battery 7.4V											
Power	110v - 240v, 50HZ / 60HZ											

AUTHORIZED
DISTRIBUTOR



WEIGHINGTEC
SYSTEMS

Call: 01823 521241
Web: weighingtec.com

Optional Accessory

Printer
Light Tower
External Display
PC Data Communicator

PRECISION BALANCE

0.01g,0.1g,1g